














Sr No	Area	Model	Product Details	Picture
1	OT	1025AX	Monitor Trolley (Economy) -Roma electrical fittings overall approx size 680mmL x700mmDx1570mmH. Suitable for monitor, camera, VCR printer, stabilizer etc. four MS shelves, MS drawer at bottom with lock. Trolley has two handles, SS IV rod, with 125mm dia Modular castors two with brakes with 08 no 3 pin plugs with powder ON/OFF switch. Complete unit powder coated finish.	
2	OT	1065XK	Bowl Stand Double two tier bowl stand approx height: 840mm with two removable SS bowl 350mm dia (304 Salem) MS tubular frame works, powder coated finish mounted on three swivel, twin wheel castors, 50mm dia.	
3	OT	1060XK	SS INSTRUMENT TROLLEY Overall approx dimensions: 29"L X 20" W X 34" H. SS tubular frame mounted on four castors, 125mm dia (Modular). Two SS brush finish shelves with three side protective railings. Attached tray is optional.	
4	OT	9001	SS Storage Cabinet for OT Equipments Approx size 36"LX18"WX60"H with two SS Shelves Complete Unit made from 304 grade Double Door with Locking systems. Floor adjustable buffer.	
5	OT	1062XK	SS DRESSING TROLLEY (LAYOUT) 1062XK Overall approx dimensions: 40"L X 22" W X 36" H. SS tubular frame mounted on four castors, 150mm dia. Two SS brush finish shelves with three side protective railings. Optional attachment SS Bowl & SS Bucket fittings	
6	OT	1044XK	SUTURES TROLLEY MODEL NO 1044XK made of complete ss sheet 304 Grade approx size 36"Lx24"Wx60"H mounted on 125 mm castors.(Trolley made of 1.2mm ss sheet).	

Sr No	Area	Model	Product Details	Picture
7	OT	OT1001	OT Instrument Trolley (OVER HEAD) made of 1.2mm stainless steel sheet & 1 ¼ ", 1.6 mm dia SS tubular frame complete unit made of 304 Salem stainless overall dimension 36"LX18"W with adjustable height from 36"to 60" with detachable SS tray size 36"x18"x2" dia with two SS side handle mounted on 125 mm modular castors two with brakes two without brake. Suitable for OT Cabin.	
8	OT	1041XK	Mayo,s Stand (1041XK) Designed for surgical procedures, sturdy and convenient, with removable stainless steel tray, adjustable height of SS tray size 560mm x400mmW, SS tubular frame mounted on four swivel twin wheel castors, 50mm dia.	
9	OT	OT1002	Hydraulic Surgeon's Chair with Back rest. Backrest adjustable made of SS tubular frame with SS round ring footrest, Cushion work on seat and back rest.	
10	OT	9003	Surgical Scrub Station with auto sensor And Knee Press overall size 72"LX30"W X 40"H made from SS tubular frame 16gauge 304 Grade & 1.6mm SS Sheet, Sink inside dimension : Top size 22", Bottom size 8" , depth Technical Specification: THREE auto shower (Jaguar fittings). Infra Ray photo sensors, Hydro Teflon Solenoids Vsize. 18" THREE auto shower (Jaguar fittings). Infra Ray photo sensors, Hydro Teflon Solenoids Hot Water Geiger 10 ltrs, Hydro polyester pump 0.5HP Hand washing Soft Shower. Electronic Circuits & Controller. GI – Pipes & Others Fittings.	 3 Bay Sink
11	OT	OT1003	SS Wash Sink Single Bay :Approx size 20"LX18"WX24"H wall mounted Complete S sheet 304 Grade	

Sr No	Area	Model	Product Details	Picture
12	OT	OT1004	SS Wash Sink Three Bay :Approx size 72"LX16"WX10"H wall mounted Complete S sheet 304 Grade72" long	
13	OT	OT1005	SS Wash Sink Single Bay :Approx size 36"LX16"WX16"H wall mounted Complete S sheet 304 Grade72" long	

Delivery within 10 to 12 weeks (Working Days) after confirm of order.

GST will be charge @ 18% extra.

Transportation To Pay Basis as actual.

Payment 60% advance with GST along with order copy

30% with GST before dispatch of the consignment 10% against installation .

Warranty Two years on Powder coating, 03 years on Welding Part, (MS tubular/CRC sheet)

02 year on Actuators (Electrical Devices)

Imported Castors and Ten Years on Stainless steel Tubular and SS Sheet of any manufacturing defects

Raw Material Details

CRCA sheet and tubular 1.6 mm 16 gauge use for Hospital Beds & Ward furniture

Powder coating process are use ASIAN PAINT pre treated through 7 tank process

Imported castors 125 mm modular use for Hospital Beds make (Tente /Setico Germany}

Stainless steel sheet are use 1.6mm (16 Gauge) 304 Grade for all SS Furniture

Stainless steel Tubular Frame are use 1 ½ x 1 ½" x1.6mm (16 Gauge) 304 Grade for all SS

Heavy duty Castors 150 mm from Kalpar Gujarat.

Fittings are also provided All SS 304 Grade in SS Furniture.

Modular Adjustable Shoe fitted on all SS Tubular frame for floor adjustment.