












| Sr No | Area    | Model                                   | Product Details   | Picture  |
|-------|---------|---|---|--|
| 1     | CAUTERY | SMART 4<br>Model 306                    | Smart 4 is a micro processor controlled based model for surgeries such as Electro vaporization Gastroenterology, cardiology, Proctology, Gynecology, Urology, Neurology, Orthopaedics & many other type of surgeries.   |    |
| 2     | CAUTERY | 400 W<br>DIGITAL<br>CAUTERY<br>MPM 333  | Tanvi MPM 333 is micro controller based model for surgeries such as Electro vapourization Gastroenterology, Cardiology, Proctology, Gynaecology, Urology, Neurology, Orthopaedics & many other types of surgeries. Digital Power adjustment for separate modes Feather touch keys are used for power setting levels in all operating mode Audio visual Alarm for disconnecting of patient plate Tissue response circuit is provided Programme selector switch is provided |    |
| 3     | CAUTERY | 300 W<br>DIGITAL<br>CAUTERY<br>MPM 3033 | We offer to our clients "Tanvi 330" units which are designed to perform both simple procedures and highly complex surgeries. Further they also provide a smooth cut, minimum tissue damage and effective coagulation.   |  |
| 4     | CAUTERY | 400 W<br>ANALOG<br>CAUTERY<br>MPM 303   | Solid State surgical Diathermy 400 w with Handswitch Pencil Analogue type   |  |
| 5     | CAUTERY | 300 W<br>ANALOG<br>CAUTERY<br>300 B     | Solid State surgical Diathermy 300 w with Handswitch Pencil Analogue type   |  |

|           |          |   |   |   |
|-----------|----------|---|---|---|
| <b>6</b>  | CAUTERY  | Radio Frequency Cautery (3MHz)          | <p>We offer to our clients superior range of <b>Radio Frequency Cautery- 3</b>, which delivers radio frequency (3MHz) wave for cutting &amp; coagulation of soft tissues. These radio waves are used for vaporizing the cells without applying any pressure on the tissue.</p>  |     |
| <b>7</b>  | OT Light | Alma Series - Tanvi 304                 | 4 reflector (heart shape)   |     |
| <b>8</b>  | OT Light | Alma Series - Tanvi 3044                | 4 + 4 reflector (heart shape)   |     |
| <b>9</b>  | OT Light | Alma Series - Tanvi 305                 | 4 reflector (Pentagonal shape)  |   |
| <b>10</b> | OT Light | Alma Series - Tanvi 3055                | 4 + 4 reflector (Pentagonal shape)  |  |
| <b>11</b> | OT TABLE | C-ARM COMPATIBLE Hydraulic Table (3015) | <p>Tanvi c-arm semi motorized table is designed in accordance with orthopedic surgeons of radiographic capability &amp; general surgery. Up and Down Movement is controlled by feather touch buttons on remote control without disturbing operating team. Radiolucent top allows clear image on c-arm machine monitor Eccentric table</p> |  |