Reagents	
Liver Function	ALT、AST、ALP、GGT、TBIL、DBIL、TP、ALB、ADA、TBA、PA、LAP、CHE
Blood Lipid	TG、CHO、LDL-C、HDL-C、ApoA-I、ApoB、Lpa、NEFA
Renal Function	UREA、CREA-E、UA、Cys-C、RBP、UMA、B2M、α-AMY
Blood Sugar	GLU-H、HbA1c、FMN
Blood Gas & Electrolyte	Ca、Co_2、Mg、Fe、P
Special Protein	ASO、RF、CRP、hs-CRP、IgA、IgG、IgM
Myocardium Zymogram	CK、CK-MB、LDH、α-HBDH、HCY



OBio BA200M
Auto Chemistry Analyzer

Technical specifications

General information	Machine type	Fully automatic analyzer, random access, STAT sample in priority
	Test speed	Chemistry: constant 200 tests/hour Chemistry+ISE(optional): constant 400 tests/hour
	Test principle	Colorimetric method, turbidimetry method, ISE method
	Test method	Endpoint method, fixed-time method (two-point method), kinetic method (rate method)
Reagent/Sample Unit	Sample/reagent tray	50~100 sample positions, support barcode scanning, automatic dilution, retesting and manual retesting. Samples can be inserted for emergency handling at any time and will be prioritized for processing. 50~100 reagent positions, support barcode scanning, 24-hour uninterrupted cooling function.
	Sample/reagent probe	1 probe, liquid level detection, volume tracking, anti-collision, clog detection(optional)
	Sample volume	Chemistry: 1.5~45μL, step by 0.1μL ISE: 70μL (serum), 140μL (urine)
	Reagent volume	10~240μL, step by 1μL
Mixing system	Mixer	1 independent mixer
Reaction unit	Cuvette	48 reaction cuvettes, each with a light diameter of 5mm
	Reaction volume	70~315μL
	Temperature	37℃±0.1℃
	Cuvette washing	2*4 stages automatic washing with alkaline cleaning solutions
Optical system	Light source	Halogen lamp, average lifespan of 2000 hours
	Wavelength	8 wavelengths, 340nm, 405nm, 450nm, 505nm, 546nm, 578nm, 630nm, 700nm
Built-in barcode scanning (optional)		Applicable to sample and reagent barcodes. Reading code formats: Codabar, ITF (Interleaved Two of Five), code128, code39, UPC/EAN, Code93
ISE module (optional)		K+,Na+,CL-
QC/Calibration	Calibration modes	K-factor method, two-point linear, multi-point linear, spline, polyline, exponential, Logi t-Log3P, Logit-Log4P, Logit-Log5P
	QC Rules	Westgard, Levey-Jennings, Cumulative sum check, Twin plot
Operation system Unit	Operation system	Windows10 64bit and above editions
	System connection	TCP/IP network connection
Others	Power Supply	100~240V, 50/60Hz
	Rated power	500VA
	Dimension	750*680*507(mm)
	Weight	65 Kg

Orchard Medical

154, Teghoria Main Road Kolkata -700157
Mobile no: **+91 84200 69980**
Website- www.orchardmedical.co.in
Email: sales.orchardmedical@gmail.com



OBio BA200M

Auto Chemistry Analyzer



OBio...
Supporting Care to Cure
A Division of Orchard Medical

OBio BA200M

Auto Chemistry Analyzer

Advanced operating system

- Optional computer with touch screen, flexible & convenient
- User friendly software interface, simple & easy to operate
- Sample automatic dilution(decrease, increase, normal)
- Reaction result automatic checking and warning, calibration result automatic checking, warning alarm log
- Linear & nonlinear calibration, multi-points calibration up to 8 points








Multifunctional Precision Sampling Probe

- Support inner and outer wall cleaning, reduce carry-over
- Liquid level detection, volume tracking, anti-collision, probe blockage alarm



Excellent Mixer Design

- Teflon coating mixer, reduce water droplet (reduce carry-over)
- Excellent mixing effect with standard mixing procedure



Reagent & Sample Tray

- 50 sample positions, it can be expanded to 100 sample positions, emergency samples in priority and can be inserted at any time
- 50 reagent positions, it can be expanded to 100 reagent positions, 24-hour continuous refrigeration function



Chemistry Reagents



Reaction Tray

- 48 recycable cuvettes
- Reaction volume 70~315 μL
- Stable & accurate temperature 37℃ for reaction



Cuvette washing

- 2*4 stages automatic washing
- Acid and alkaline cleaning solutions and has a residual alarm function



OBio Intelligent-Sampling technology creates an innovative platform for HbA1c detection, supporting whole blood testing, easy to use and highly accurate.

Convenient operation

No pretreatment is required, the whole blood sample can be directy loaded onto the machine
Intelligent-Sampling technology supports the absorption of RBC below the sample liquid surface

Accurate and reliable

The results are traceable. The Intelligent-Sampling sampling technology avoids errors caused by insufficient effective RBC sampling volume and provides reliable diagnostic data for clinical practice

